

RBLLRRNoneU41Y
(Book 33)

M4232

M4233

M4234

M4235

PS-16738/KLH hENOS ϵ -NH₂ terminal
ASPWKFLFKTKFHQYAMNGE/KLH conjugate

- ① 1st immunization I.P. ϵ FCA @ 50 μ g/mouse 11/18/95 RW
- ② 2nd immunization I.P. ϵ FCA @ 50 μ g/mouse 12/7/95
- ③ TEST BLEED 12/28/95
- ④ TEST BLEED 03/04/96
- ⑤ KK mice 3/29/96

Test Bleed of M #4232-35

EASY BEAM

Plate # 2 Date 12/30/85 Filter 492 nm
 Operator RW Comment PS-1673B @ 100ng/well hiNOS @ -NH₂ Term

blank 1 100 2 200 3 400 4 800 5 1000 6 3200 7 6400 8 12800 9 25600 10 51200 11 102400 12

M4232

A	+0.000	+0.000	+0.000	+0.005	+0.005	+0.000	+0.005	+0.000	+0.000	+0.000	+0.000	+0.000	A
---	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	---

33

B	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	B
---	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	---

34

C	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	C
---	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	---

35

D	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.000	+0.018	+0.000	+0.000	+0.000	D
---	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	---

M4232

E	+0.000	+0.021	+0.008	+0.000	+0.000	+0.000	+0.000	+0.000	+0.006	+0.015	+0.000	+0.000	E
---	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	---

33

F	+0.000	+0.023	+0.017	+0.010	+0.000	+0.000	+0.013	+0.000	+0.000	+0.000	+0.000	+0.000	F
---	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	---

34

G	+0.000	+0.039	+0.023	+0.011	+0.008	+0.000	+0.006	+0.011	+0.000	+0.000	+0.000	+0.000	G
---	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	---

35

H	+0.000	+0.159	+0.036	+0.015	+0.010	+0.000	+0.010	+0.000	+0.000	+0.000	+0.000	+0.000	H
---	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	---

1 2 3 4 5 6 7 8 9 10 11 12

M4232 - 4235 Test Bleed

EASY BEAM

Plate # 2 Date 3/5/96 Filter 492 nm
 Operator RW Comment PS-1673B @ 100ng/well HIVOS NH₂ term

Blank 1 100 200 400 800 1600 3200 6400 12800 25600 51200 102400

M4232
 A +0.000 +0.039 +0.012 +0.000 -0.006 -0.008 -0.005 +0.006 +0.000 +0.015 +0.000 +0.000 A

4233
 B +0.000 +0.009 +0.000 +0.000 +0.000 +0.000 +0.006 +0.006 +0.000 +0.000 +0.000 +0.000 B

4234
 C +0.000 +0.036 +0.014 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 C

4235
 D +0.000 +0.025 +0.005 +0.031 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 D

M4232
 E +0.000 +0.017 +0.008 +0.000 +0.000 +0.000 +0.000 +0.000 +0.019 +0.000 +0.000 +0.000 E

4233
 F +0.000 +0.024 +0.008 +0.000 +0.000 +0.000 +0.008 -0.005 +0.000 +0.000 +0.000 +0.000 F

4234
 G +0.000 +0.009 +0.000 +0.000 +0.000 +0.000 +0.011 +0.000 +0.000 +0.000 +0.000 +0.000 G

4235
 H +0.000 +0.016 +0.008 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 +0.000 H

1 2 3 4 5 6 7 8 9 10 11 12

No Ab :. KK

BEST AVAILABLE COPY